

Oxygen & Suction Therapy

Vacuum Regulator

Silbermann's vacuum regulator is intended to measure and adjust the vacuum level generated from the central vacuum system.

Its robust design insures a safe, reliable and continuous operation that is planned to last.

The vacuum regulator includes an inlet male adaptor, shutoff valve, 0-760 or 0-200 mmHg gauge and vacuum DISS male thread inlet connection.





Features & Benefits

- Standard model vacuum scale:0-760mmHg
- Low Vacuum model vacuum scale:0-250 mmHg
- Both types with Vacuum DISS male outlet and 1/8" NPT inlet female
- On/Off switch valve
- Rail mounted regulator includes
 1 meter flexible medical tubing and
 probe, and a universal rail clamp

Silbermann's Quality Management System:

- EN ISO 13485-2016
- ISO 9001-2015
- CE marked







Catalogue Numbers - Vacuum Regulator with probe



Type of Male Adapter	0-760 mmHg	0-250 mmHg
Silbermann- Ohio	4602000	4606000
ENV 737-6	4602100	4605000
DIN 13260-2	4668000	4669000
Puritan Bennett type	4604000	4608000
BS5682	4602150	**TBD
Diss	4611000	4611500
1/8 NPT female thread	4012000	4013000



Catalogue Numbers - Vacuum Regulator on a hanging rail with probe*



Type of Male Adapter	0-760 mmHg	0-250 mmHg
Silbermann- Ohio	4621000	4625000
ENV 737-6	4621100	4621200
DIN 13260-2	4672000	4637000
Puritan Bennett type	4623000	4627000
BS5682	**TBD	**TBD
Diss	4620000	4624000

^{*} Includes 1 meter flexible medical tubing and probe, and a universal clamp.

^{**} TBD -Available on customer request





Assembly - Suction Unit

The suction unit consists of:

- Vacuum Regulator with probe
- Polysolfone overflow safety trap bottle (P.N 4034000)
- 2 liter polycarbonate collecting bottle (P.N 4053000)*
- Collecting jar stainless steel holder (P.N 4000253)
- 1 meter P.V.C Vacuum hose (P.N 2118100)
- 2 meter P.V.C Vacuum hose (P.N 2118200)
- Hydrophobic filter (P.N 4415312) Optional



Suction Unit on Hanging Rail

Suction unit can be supplied rail mounted or direct mounted.

* Optional - 2 liter polysolphone collecting bottle



Related Products



Vacuum slide with hanging rail clamp P.N - 3076100



Vacuum slide for wall mounting or BHU mounting P.N - 3076000

Distribu	ited B	у:		

SILBERMANN Technologies Ltd.

Address: #5 Harakevet st. POB 4605, Petah-Tikva, 4900715, ISRAEL

Tel: +972-3-9309186 | Fax: +972-3-9315721 E-Mail: info@silbermnn.com, www.silbermnn.com

